



京天成生物技术（北京）有限公司

AbMax Biotechnology Co., LTD

北京海淀区上地东路 29 号 2B 702

Building-2B, Suite 702, 29 E ShangDi Road,

Beijing, 100085, China

Tel: 86-10-82896224, Fax: 86-10-62977292

E-mail: info@antibodychina.com

www.antibodychina.com

Anti-Rabies Virus, Clone 5A2

(mouse monoclonal IgG), FITC-conjugated

Catalog # 12-0018

Lot # 100021

Antibody Class: IgG mouse monoclonal antibody produced *in vitro* by mouse-mouse hybridoma 5A2, purified by Protein G-Sepharose chromatography and cross-linked to FITC.

Immunogen: Rabies virus

Formulation: 200 reactions. IgG-FITC conjugate in 500µl PBS.

Physical State: Frozen liquid

Storage: Store the product away from light at -20°C. Product is stable for about 6 months at 2-8°C. For maximum recovery of product, centrifuge the original vial prior to removing the cap.

FOR IN VITRO RESEARCH USE ONLY

NOT FOR USE IN HUMANS

Research Applications

Immunohistochemistry: use at 1:20 dilution.
