



京天成生物技术（北京）有限公司

AbMax Biotechnology Co., LTD

北京海淀区上地东路 29 号 2B 702

Building-2B, Suite 702, 29 E ShangDi Road,

Beijing, 100085, China

Tel: 86-10-82896224, Fax: 86-10-62977292

E-mail: info@antibodychina.com

www.antibodychina.com

Anti-HIV/p24, Clone 1G12

(mouse monoclonal IgG)

Catalog # 05-0013

Lot # 100077

Immunogen: Synthetic peptide corresponding to the middle region of p24 protein of HIV.

Specificity: Recognizes p24 of HIV.

Species Cross-reactivity: Not determined.

Formulation: 0.2mg of mouse monoclonal IgG in 0.2ml of 0.01M Tris-HCl, pH 8.0, 0.15M NaCl, 0.02% sodium azide. Protein was determined by OD₂₈₀ absorbance.

Physical State: Frozen liquid

Storage: Store the product at -20 °C . Product is stable for about 6 weeks at 2-8 °C as an undiluted liquid. Prepare working dilution fresh each day. Avoid repeated freezing and thawing.

**FOR IN VITRO RESEARCH USE ONLY
NOT FOR USE IN HUMANS**

Research Applications

Sandwich ELISA: Used as capturing antibody, 5~10µg/ml.

Indirect ELISA: 0.5~2µg/ml. Detection sensitivity 3ng/ml.

Immunoblot Analysis: 0.5~2µg/ml.

