



京天成生物技术（北京）有限公司

AbMax Biotechnology Co., LTD

北京海淀区上地东路 29 号 2B 702

Building-2B, Suite 702, 29 E ShangDi Road,

Beijing, 100085, China

Tel: 86-10-82896224, Fax: 86-10-62977292

E-mail: info@antibodychina.com

www.antibodychina.com

## Anti-HBsAg, Clone 6D9

(mouse monoclonal IgG)

Catalog # 05-0022

Lot # 100016

**Immunogen:** Recombinant human Hepatitis B Surface Antigen (HBsAg), Adw

**Specificity:** Recognizes human HBsAg(Adw and ADr)

**Species Cross-reactivity:** Not determined

**Formulation:** 0.2mg of mouse monoclonal IgG in 0.2ml of 0.01M Tris-HCl, pH 8.0, 0.15M NaCl, 0.02% sodium azide. Protein was determined by OD280 absorbance.

**Physical State:** Frozen liquid

**Storage:** Store the product at  $-20^{\circ}\text{C}$ . Product is stable for about 6 weeks at  $2-8^{\circ}\text{C}$  as an undiluted liquid. Prepare working dilution fresh each day. Avoid repeated freezing and thawing.

**FOR IN VITRO RESEARCH USE ONLY  
NOT FOR USE IN HUMANS**

---

### Research Applications

**Indirect ELISA:** A working concentration of 0.5~2 $\mu\text{g/ml}$  is suggested. Detection sensitivity 0.1ng/ml.

**Sandwich ELISA:** Used as capturing antibody, 5~10 $\mu\text{g/ml}$ .

---