



京天成生物技术（北京）有限公司

AbMax Biotechnology Co., LTD

北京海淀区上地东路 29 号 2B 702

Building-2B, Suite 702, 29 E ShangDi Road,

Beijing, 100085, China

Tel: 86-10-82896224, Fax: 86-10-62977292

E-mail: info@antibodychina.com

www.antibodychina.com

Anti-*E.coli* M15 HCP Antibody

(rabbit polyclonal IgG), HRP-conjugated

Catalog # 12-0005

Lot # 100091

Antibody Class: Rabbit polyclonal antibody IgG, purified by Protein G-Sepharose chromatography and cross-linked to horseradish peroxidase.

Immunogen: *E.coli* M15 HCP(host cell protein)

Formulation: 0.2mg of IgG-HRP conjugate in 0.2ml of PBS, 25% glycerol.

Physical State: Frozen liquid

Storage: Store the product away from light at -20°C.

Product is stable for about 6 months at 2-8°C. For maximum recovery of product, centrifuge the original vial prior to removing the cap.

FOR IN VITRO RESEARCH USE ONLY

NOT FOR USE IN HUMANS

Research Applications

Direct ELISA or Sandwich ELISA:

Used as detecting antibody, a working dilution of 1:500 to 1:1000 is suggested.
