



京天成生物技术（北京）有限公司

AbMax Biotechnology Co., LTD

北京海淀区上地东路 29 号 2B 702

Building-2B, Suite 702, 29 E ShangDi Road,

Beijing, 100085, China

Tel: 86-10-82896224, Fax: 86-10-62977292

E-mail: info@antibodychina.com

www.antibodychina.com

Anti-HIV/p24, Clone 1D4
(mouse monoclonal IgG) , HRP-conjugated

Catalog # 12-0011

Lot # 100084

Antibody Class: mouse monoclonal antibody IgG produced *in vitro* by mouse-mouse hybridoma 1D4, purified by Protein G-Sepharose chromatography and cross-linked to horseradish peroxidase.

Immunogen: Synthetic peptide corresponding to the NT of p24 protein of HIV.

Formulation: 0.2mg of IgG-HRP conjugate in 0.2ml of PBS, 25% glycerol.

Physical State: Frozen liquid

Storage: Store the product away from light at -20°C. Product is stable for about 6 months at 2-8°C. For maximum recovery of product, centrifuge the original vial prior to removing the cap.

FOR IN VITRO RESEARCH USE ONLY
NOT FOR USE IN HUMANS

Research Applications

Direct ELISA or Sandwich ELISA:

Used as detecting antibody, a working dilution of 1:500 to 1:1000 is suggested.
