



京天成生物技术（北京）有限公司

AbMax Biotechnology Co., LTD

北京海淀区上地东路 29 号 2B 702

Building-2B, Suite 702, 29 E ShangDi Road,

Beijing, 100085, China

Tel: 86-10-82896224, Fax: 86-10-62977292

E-mail: info@antibodychina.com

www.antibodychina.com

Anti-Rabies Virus, Clone 1C7

(mouse monoclonal IgG)

Catalog # 05-0020

Lot # 100020

Immunogen: Rabies virus

Specificity: Recognizes rabies virus.

Formulation: 0.2mg of mouse monoclonal IgG in 0.2ml of 0.01M Tris-HCl, pH 8.0, 0.15M NaCl, 0.02% sodium azide. Protein was determined by OD280 absorbance.

Physical State: Frozen liquid

Storage: Store the product at -20°C . Product is stable for about 6 weeks at $2-8^{\circ}\text{C}$ as an undiluted liquid. Prepare working dilution fresh each day. Avoid repeated freezing and thawing.

**FOR IN VITRO RESEARCH USE ONLY
NOT FOR USE IN HUMANS**

Research Applications

Indirect ELISA: A working concentration of 0.5~2 $\mu\text{g}/\text{ml}$ is suggested.

Immunohistochemistry: A working concentration of 2~10 $\mu\text{g}/\text{ml}$ is suggested

Sandwich ELISA: Used as capturing antibody. 5~10 $\mu\text{g}/\text{ml}$.
