



京天成生物技术（北京）有限公司

AbMax Biotechnology Co., LTD

北京海淀区上地东路 29 号 2B 702

Building-2B, Suite 702, 29 E ShangDi Road,

Beijing, 100085, China

Tel: 86-10-82896224, Fax: 86-10-62977292

E-mail: info@antibodychina.com

www.antibodychina.com

## Anti-Human Troponin T Antibody

(rabbit polyclonal IgG)

Catalog # 06-0006

Lot # 100082

**Immunogen:** Recombinant Human Troponin T

**Specificity:** Recognizes human Troponin T. Troponin T has been shown to be a Biomarker for heart damage.

**Species Cross-reactivity:** Not determined.

**Formulation:** 0.2mg of rabbit polyclonal IgG in 0.2ml of 0.01M Tris-HCl, pH 8.0, 0.15M NaCl, 0.02% sodium azide. Protein was determined by OD<sub>280</sub> absorbance.

**Physical State:** Frozen liquid

**Storage:** Store the product at -20 °C . Product is stable for about 6 weeks at 2-8 °C as an undiluted liquid. Prepare working dilution fresh each day. Avoid repeated freezing and thawing.

**FOR IN VITRO RESEARCH USE ONLY  
NOT FOR USE IN HUMANS**

---

### Research Applications

**Indirect ELISA:** A working concentration of 0.1~2µg/ml is suggested.

**Sandwich ELISA:** Used as capturing antibody, a working concentration of 5~10µg/ml is suggested.

**Immunoblot Analysis:** 0.5~2µg/ml.

---