



京天成生物技术（北京）有限公司

AbMax Biotechnology Co., LTD

北京海淀区上地东路 29 号 2B 702

Building-2B, Suite 702, 29 E ShangDi Road,

Beijing, 100085, China

Tel: 86-10-82896224, Fax: 86-10-62977292

E-mail: info@antibodychina.com

www.antibodychina.com

## Anti-Human IgG Antibody

(rabbit polyclonal IgG), HRP-conjugated

Catalog # 12-0030

Lot # 100074

**Antibody Class:** Rabbit polyclonal antibody IgG, purified by Protein G-Sepharose chromatography and cross-linked to horseradish peroxidase.

**Immunogen:** Human IgG

**Formulation:** 0.2mg of IgG-HRP conjugate in 0.2ml of PBS, 25% glycerol.

**Physical State:** Frozen liquid

**Storage:** Store the product away from light at -20°C.

Product is stable for about 6 months at 2-8°C. For maximum recovery of product, centrifuge the original vial prior to removing the cap.

**FOR IN VITRO RESEARCH USE ONLY**

**NOT FOR USE IN HUMANS**

---

### Research Applications

**Direct ELISA or Sandwich ELISA:**

Used as detecting antibody, a working dilution of 1:1000 to 1:2000 is suggested.